
amtech®

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Systems, facilities and equipment  
for cleanrooms, assistance, calibration  
and validation services

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AM Instruments offers systems, facilities and equipment for cleanrooms. AMTech® manufactures standard and custom-made products through a careful analysis of the customer's requests and special needs. Our team is able to offer solutions even for complex dimensions, functions and layouts thanks to the support of the highly modern design and production processes.

FROM DESIGN TO PROJECT TO REALISATION, AM INSTRUMENTS' APPROACH IS PROFESSIONAL, COMPETENT AND PROACTIVE.

With its own multifunctional laboratory featuring dedicated working areas for the maintenance, repair and calibration of instruments, AM Instruments offers assistance and services which are among the most professional on the market.

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## SERVIZI

- ASSISTENZA E TARATURA
- CONVALIDA

## SYSTEMS AND FACILITIES FOR CLEANROOMS

- LAF SYSTEMS
- FACILITY MONITORING SYSTEMS
- EQUIPMENT

## SERVICES

- TECHNICAL ASSISTANCE & CALIBRATION
- VALIDATION

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
AM Instruments has designed the standard **Multiplo®** module, a stand-alone laminar flow unit available in four sizes and assembled according to the customer's specific needs.

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**AM Instruments** designs and builds LAF systems for contamination controlled areas, in accordance with applicable legislation (ISO 14644-1).

A team of engineers, technicians and qualified operators support the customer from development to delivery, and provide after-sales assistance.

Advanced technology and design, adaptability to spaces and operational excellence.

Sectors: pharmaceutical, chemical, electronic and food; in production and research processes.

- laminar flow modules, assembled to suit different dimensions and layouts
- down-cross units
- traditional and new generation pass-boxes with HPE technology
- RABS
- battery-powered LAF trolleys designed for the transfer of containers holding potentially exposed product between Grade A areas, passing through areas with a lower classification in full safety
- air showers
- LAF cabinets

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The Multiplo® module is a self-contained LAF unit with a HEPA filter and fan, which comes in four standard sizes. The machine can be assembled to suit a variety of sizes and layouts to meet class-ISO5 requirements (compliant with ISO 14644-1) of certain sectors and specific phases of production. The standard module is made entirely from AISI 304 stainless steel with a Scotch Brite brushed finishing, and is fitted with:

- EU14 HEPA filter - 99.995 MPPS efficiency (particles 0.12 - 0.24  $\mu\text{m}$ ),  $\geq$  99.999 DOP (0.3  $\mu\text{m}$ )
- centrifugal motor fan, IP55 protection rating
- plenum made of stainless steel, designed for 100% DOP entry

Modules can be supplemented with a wide range of accessories.

## WORKING PRINCIPLE

Air is drawn from the environment through pre-filters located at the sides then forced by the motor fan into the plenum, with optimised structural sections, which distributes pressure evenly on the surface of the HEPA filter in order to facilitate laminar airflow downstream.

The areas surrounding the motor fan are kept at negative pressure, thus creating a dynamic seal and enhancing the intrinsic safety of the Multiplo® system.

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Design and production of entirely custom-made LAF modules, with customised solutions based on individual applications. Self-regulating electronic motors guarantee constant air speed for the entire lifespan of the absolute filter, without the intervention of the operator. The modules guarantee reduced maintenance and ease of integration in centralised control systems.

**Maximised laminar flow areas:** with the new equaliser cover, AM Instruments eliminates the non-fluxed areas and edges between HEPA filters, reduces the use of special filters and makes it possible to make generate flux even in areas that are not perfectly rectangular. Filtering coverage per single module is the highest that has ever been achieved. The integration of LED bar lights inside the module makes all surfaces smooth, eliminating areas that are difficult to clean.

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The AM Instruments down-cross booth is a vertical laminar flow unit for weighing or sampling rooms. It guarantees operator exposure levels of  $<100\mu\text{g}/\text{m}^3$  TWA, which can be further reduced with the use of flexible or rigid barriers to  $<10\mu\text{g}/\text{m}^3$  TWA.

- top quality finishing for the pharmaceutical sector
- cGMP compliance
- high-performance fans
- AISI 304 stainless steel structure
- integrated electrical panel
- alarm management
- validation
- IQ and OQ documentation in accordance with cGMP standards
- custom-made solutions

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## A NEW GENERATION OF FLUXED BIODECONTAMINATION PASS-BOXES

AM Instruments has applied the most innovative technologies to develop a new pass-box model.

- fluxed pass-box and pass-boxes, stand-alone or integrated into existing ventilation systems
- material transfer chambers within adaptable shapes that can be adjusted to suit every need
- solutions also proposed with HMI touch screen operator panels for the traceability of transfer cycles and door opening times
- possibility of option to integrate UV lamps and hydrogen peroxide decontamination systems.

Using integrated low-concentration hydrogen peroxide systems (<8%) with HPE technology, our decontamination chambers guarantee a degree of bacterial reduction equal to log 6 in extremely short times. Ideal for the decontamination of thermolabile materials and electronic components.

MyBox® is produced in accordance with the GMP standards for pharmaceutical environments, with high-quality and high-finish materials. It can perform automatic biodecontamination cycles through an integrated PLC, adapting to the specifications of the individual load.

During the validation phase, chemical and microbiological mapping with bioindicators (*Geobacillus Stearothermophilus*) ensures the correct distribution and effectiveness of the hydrogen peroxide in guaranteeing the desired level of biodecontamination. MyBox® can be easily installed for the passage of material from classified C/D areas into sterile Grade B areas, guaranteeing the sterilisation of external surfaces. This device is essential when the products/components to be sterilised are not autoclavable.

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## GENERATOR WITH HPE TECHNOLOGY

HPE technology is based on a hydrogen peroxide solution with a concentration lower than 8% which, passing through a high voltage electromagnetic field, undergoes an ionisation process allowing a high and effective diffusion with validated biodecontamination cycles and a reduction of the surface microbial load of log 6, including viruses, bacteria, moulds and microbial spores. The system is equipped with a control panel that can be connected to a supervision system for complete automation of the cycle.

In the environment, the ionised mist behaves with the dispersion characteristics of a gas.

- highly compatible with all materials used in controlled environments
- utilities required: compressed air, connections to an electrical supply and to an external exhaust system
- simple and comprehensive design suitable for any need

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The drive towards the implementation of more modern, safer technologies in aseptic pharmaceutical production encourages production sites and processes to change rapidly to adapt to the increasingly stringent demands of regulatory bodies. The use of physical separations (Restricted Access Barrier Systems) between the operator and the product must be considered a primary solution in the design of new filling lines or the renovation of existing ones.

AMTech® offers feasibility studies to align production lines with the most innovative design solutions. There is the option to design custom-made RABS with glove doors positioned in accordance with the specific needs of each process.

Transfer ports are used to eliminate contamination due to the transfer of critical process materials through Grade B areas into Grade A areas.

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## CLASS ISO 5 MATERIAL TRANSFER SYSTEM

The LAF trolley has been developed to allow the completely safe handling and transfer of materials between Grade A areas, passing through a Grade B area, without the need for an electrical supply (battery life of up to one hour).

The LAF trolley is commonly used to transfer partially stoppered vials from the filling line to the freeze-dryers, and then from the freeze-dryers to the vial-capping machines. Air is drawn from the environment through pre-filters located at the sides, filtered by the HEPA filter, forced horizontally into the working area then expelled through openings located on the doors.

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## MATERIAL TRANSFER SYSTEMS WITH CONTINUOUS PARTICLE MONITORING


New design with continuous particle counter:

- sealed doors and recirculating airflow
- sampling probe in the "worst case" position close to recovery
- compliance with data integrity in the management of monitoring data

- complete control of process parameters through the integrated comfort panel

"No continuous monitoring of particles is carried out throughout the critical process for A Grade areas. Transportation of not finally stoppered vials with products through B Grade area is carried out in portable untight carts with horizontal Grade A LAF, no continuous particle monitoring is carried out in the carts. The particles are checked only before transportation beginning and after product unloading."

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## AIRFLOW DECONTAMINATION SYSTEM

System for the decontamination of personnel entering a contamination-controlled area. The shower is designed for the entrance and exit of one person at a time and is equipped with two access doors, normally open, which are interlocked so as not to allow simultaneous opening. A light signal on each side displays the status of the doors. Entering the air shower from the "unclassified" side automatically activates the airflow when the door closes, according to a pre-set operating cycle.

Emergency buttons located inside the air shower and externally on both sides allow the system to be halted immediately, opening the doors. The air shower can also be used for materials on pallets being transferred from an unclassified area to a classified one (e.g. Grade D). The interior is kept at constant negative pressure in relation to the adjoining rooms by means of a fan that expels air.

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The laminar flow cabinets are made of AISI 304 stainless steel and are equipped with ventilation and air filtering systems to guarantee ISO 5 conditions inside the cabinet when the doors are closed. Scotch Brite finishing (polished upon request) for both external surfaces and the internal compartment. Fully customisable sizes.

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AM Instruments designs and produces particle and microbial monitoring systems with integrated automation, which comply with cGMP, GAMP and ISO guidelines.

A single contact person will handle all phases of the process, from initial inspections to the planning and design of the system, from the creation to integrated software developed according to the company's needs to installation and inspections and from validation to operator training.

### COMPLETELY CUSTOM-MADE


- feasibility analysis in the pre-sale phase
- a design and development team entirely at your service
- no limit to integration with existing machines and software

### IMMEDIATE

- immediate reporting
- complex reports obtainable with a single click
- macro-operations translated into simple software operations

### SECURE

- compliance with data integrity
- system controlled accesses
- reporting compliant with Federal Register 21CFR Part 11 and EudraLex Vol. 4 Annex 11 regulations on electronic signature management and electronic records.
- system-controlled accesses

### COMPLIANT

- In line with the following standards:
  - > cGMP,
  - > GAMP 5,
  - > 21 CFR Part 11,
  - > Eudralex Annex 1, Annex 11,
  - > ISO 14644-1

### REAL-TIME ASSISTANCE

- robust and reliable systems
- ongoing assistance – including remote assistance – from our technicians, with real time interventions
- predictive maintenance, which allows system performance to be measured and maintenance work to be carried out only when it is actually required

### PROJECT PLAN AND VALIDATION

The design and implementation of continuous monitoring systems are accompanied by validation documentation, as required by current regulations.

- draft of the system's functional specifications
- project plan
- program specifications
- hardware specification/system installation details
- structural notes
- system testing /FAT
- hardware and software design
- IQ/OQ
- user manuals
- staff training certificates
- change control/error reports
- validation report

## COMPONENTS

- remote particle counters or counters with incorporated pumps
- sampling devices with SmartProbe® automatic closing system (AM Instruments patent)
- membrane or impact microbial sampling heads
- electrical control panel
- SCADA supervision system
- interfacing with process and production machines

### AFTER SALES ASSISTANCE

|                                   |   |
|-----------------------------------|---|
| Periodic calibrations             | Interventions on request                                    |
| Planned maintenance interventions | Calibration on customer's premises and/or at AM Instruments |
| Tailored maintenance contracts    | Telephone and/or remote consulting and assistance           |

### TYPES OF CONTINUOUS PARTICLE MONITORING SYSTEMS

|  |  |
|--|--|
| Single-sensor systems with central aspiration unit | Single-sensor systems with built-in pump |
|--|--|


### MICROBIAL CONTROL SYSTEM MEASUREMENT TECHNOLOGIES

|                       |
|-----------------------|
| Gel membranes         |
| Impact sampling heads |

### ENVIRONMENTAL PARAMETERS DETECTABLE BY THE SYSTEMS

|             |           |
|-------------|-----------|
| Temperature | Humidity  |
| Pressure    | Air speed |

All analogue, digital, serial, and/or Ethernet signals generated by the section and/or process (e.g. sanitisation stage active, HVAC anomaly, alarm for the presence of solvents)

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## FOG GENERATOR FOR THE VISUAL CONTROL OF AIR FLOWS IN CLASSIFIED ENVIRONMENTS

Thanks to the standard remote control option, a single operator can switch the unit on/off and adjust the intensity and speed of the mist. Works with room temperature process water. MyFog® has an intuitive interface with a 2.8" touch screen that allows the operator immediate access to visualisation and diagnostics of the main operating parameters.

- visualisation of the speed and direction of air flows in all classified environments
- balancing of pressures between rooms via visual indicator
- support in detecting air stagnation areas

## TECHNICAL SPECIFICATIONS

- remote control via RC (fan speed, fog density, pause mode)
- 2.8" Tft touch screen
- water level visualisation and diagnostics with colour indication of status
- temperature monitoring and diagnostics with graphic visualisation
- automatic optimisation of the life cycle of transducers, which are switched on selectively depending on the hours of operation
- fog density adjustment
- fan rotation speed adjustment (fog diffusion speed)
- fast Start/Stop function
- monitoring of operation times, emission, and use of every single piezoelectric transducer
- piezoelectric transducer malfunction diagnostics
- explicit alarm diagnostics indicator
- diagnostics menu in operator mode (read-only) with password-protection
- diagnostics menu in supervisor mode (changes allowed) with password-protection
- pause function with storage of current generation settings
- cycle resumption with stored settings upon exit from pause or after power cut
- variable acoustic signals depending on the function
- transducer lifecycle alarm

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## BROADCAST MESSAGING SYSTEM FOR CLEANROOMS

MyPage+®, which is designed specifically for cleanroom environments, offers the option to display a variety of templates and a number of different messages at the same time. Messages can also be scheduled and a time limit can be set after which a message will disappear.

There is no limit to the number of devices that can be installed inside a facility. MyPage+® can be managed individually or in groups, assigning a specific ID to each of them, and allowing messages to be sent simultaneously to the various areas of the plant.

### FEATURES AND BENEFITS

- displays scheduled or instant messages
- displays times based on the NTP (network time protocol)
- displays one or more templates containing videos, images, text
- displays emergency messages with priority over other content
- option to add content via web application with any electronic device logged on to the network
- access to the application via username and password, with different levels of authorisation
- option to manage MyPage+® as an individual device or in groups of devices
- option to connect to the customer's SCADA system to show values such as differential room pressure, temperature, humidity, particle count etc.

# ASSISTENZA E TARATURA TECHNICAL ASSISTANCE AND CALIBRATION

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AM Labo's multifunctional laboratory combines the best technology with our technicians' skills and experience. Calibrations, repairs and maintenance of the instruments entrusted to us are performed in controlled-contamination areas.

- spare parts warehouse
- warehouse for incoming/outgoing instruments
- inlet air filtration: HEPA H14
- control of temperature and relative humidity
- anti-static floor and ESD work areas

With its own multifunctional laboratory featuring dedicated working areas for the calibration of instruments and research and development support activities, **AM Instruments** offers assistance and services which are among the most professional on the market.

- air and liquid particle counters
- Lighthouse, HACH, Pamas, Royco, MetOne, Hiac

## SERVICES

- after-sales assistance for ordinary and extraordinary maintenance of the instruments

- calibration of particle counters in accordance with ISO 21501-4 and ISO 21501-3 standards
- on-site calibration: continuous monitoring instruments and systems are calibrated at the customer's site
- wait&go: calibration of the instruments in real time with a simple appointment
- upon request a replacement counter provided for the length of time needed for repair or calibration

## DOCUMENTATION

- calibration report: the instrument is calibrated in accordance with the manufacturer's procedures, using particles of certified dimensions as defined by ISO 21501-4, ISO 21501-3 and JIS 9921 B-1997 and instruments pursuant to LAT or equivalent standards
- technical assistance report: description of the activities carried out, including any suggestions and training for the correct use of the instruments

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# VALIDATION

---

The group of technicians and specialists that make up the AM Instruments validation team is able to offer a wide range of services that include not only the validations of HVAC and LAF systems but also those of autoclaves, ovens, depyrogenating tunnels, dryers and process utilities. The team's activities include calibration and operation verification of critical instrumentation for demi, PW and WFI water generation plants, risk analysis.

## DOCUMENTATION

The drafting of documents and DQ, IQ, OQ, PQ, FAT, SAT protocols, in compliance with sector regulations, is performed by experts whose task is to produce documentation in line with the requests of national and international inspection authorities.

## TRAINING AND ADVICE

The AMTech® validation team organises courses on request at your site to train the staff involved in the use and maintenance of the tools, with targeted insights into your problems.

Through our consulting service, we can offer our support in the choice of measuring instruments, evaluating the level of accuracy necessary for your processing, the best value for money and optimisation of the available machinery.

---

### HVAC LAMINAR FLOWS

- drafting and execution of protocols
- checking the speed, flow rate and uniformity of the airflow
- checking the integrity of HEPA or ULPA filters
- checking the differential pressures between the premises
- checking the flow of air (Smoke Test)
- checking the thermohygrometric conditions
- checking the environmental classification
- checking recovery time

### EQUIPMENT UTILITIES

- drafting and execution of protocols
- autoclaves
- incubators
- climatic chambers
- ultrafreezers
- refrigerators
- depyrogenation ovens
- de-aspirator tunnel
- PW, WFI
- nitrogen
- compressed air
- steam quality
- clean steam generator
- cryogenic refrigerators

### INSTRUMENTS CALIBRATION

- pressure gauges
- temperature probes
- timers
- CO<sub>2</sub>
- speed

## PROCESS CONTROL

**CALIBRATION, MAINTENANCE AND FUNCTIONAL CHECK ON ANY MEASURING INSTRUMENT OF PHYSICAL AND CHEMICAL PARAMETERS. MEASURING YOUR PROCESS RELIABLY AND ACCURATELY ALLOWS YOU TO WORK AND PRODUCE IN SAFETY AND QUALITY, ALSO OBTAINING A SIGNIFICANT REDUCTION IN COSTS.**

### PHYSICAL PARAMETERS

#### TEMPERATURE

- PT100, PT1000, thermocouples, thermometers, thermostats
- working range from -50 to 1200 ° C (please contact the team for other needs specifications)

#### PRESSURE

- absolute, relative and differential, vacuum gauges, pressure switches
- working range from -980 mBar to 100 Bar (please contact the team for other needs specifications)

#### LIQUID FLOW

- volumetric flow rate in full piping, flow rate in open channel
- working range from DN08 to DN300 for flow rate in full piping, no limitation for flow rate in open channel

#### LEVEL

- ultrasound, radar, capacitive rod, hydrostatic head

### CHEMICAL / PHYSICAL PARAMETERS

- pH
  - redox
  - conductivity: from concentration / high field measurements to the measures for pure / ultrapure water according to USP / EP standard
  - dissolved oxygen
  - turbidity
  - total free chlorine
  - hardness
  - additional parameters: ammonia, nitrates, phosphates, TOC, COD, chromium, iron, manganese, aluminum, copper and silicates
- OTHER INSTRUMENTS: AUTOMATIC SYSTEMS FOR WATER SAMPLING.

## APPLICATIONS

- pharmaceutical
- food and beverage
- chemistry
- chemical-pharmaceutical
- hydraulics

f

in

w

COMPANY WITH  
QUALITY SYSTEM  
CERTIFIED BY DNV GL  
= ISO 9001:2015 =

03-024-21 rev 02

**aminstruments**

AM Instruments srl  
Via Isonzo, 1/C  
20812 Limbiate (MB) - Italy  
Tel. +39 02 872892.1

[aminstruments.com](http://aminstruments.com)  
[info@aminstruments.com](mailto:info@aminstruments.com)