

BIOCLEAN FUSION Sterile Polychloroprene Gloves S-BFAP

About

BioClean Fusion™ Sterile Polychloroprene Cleanroom Gloves contain no natural latex proteins, are comfortable and thin to reduce wearer fatigue and increase tactility. The antistatic glove is doubledonnable and features a beaded cuff for strength.

Specifications

COMPATIBILITY:

ISO Class 4 & EU GMP Grade A

LENGTH:

290mm (11.4")

MATERIAL:

Polychloroprene

PROTEIN LEVEL:

Polychloroprene contains no natural latex proteins

SURFACE: Finger-textured

SHAPE:

Ambidextrous

COLOUR:

Green

STERILIZATION:

Gamma irradiation

Features

- Increased tactility
- Powder-free & latex-free
- Antistatic
- · Beaded cuff for strength
- Non-particulating EasyTear[™] packaging







Quality Standards

- Manufactured in a facility holding ISO 9001: 2008 and 89/686/EEC
 Article 11b certifications
- Conforms to Category 3 Complex Design Personal Protective Equipment PPE Directive 89/686/EEC
- Processed in an NEBB certified ISO Class 4 environment
- Complies with the requirements of Glove Standard EN420:2003 and EN374-1, 2 & 3:2003
- Physical properties comply with European Medical Glove Standard EN455-2:2009
- Sterilized in accordance with ISO 11137-1:2006
- Suitable for contact with all foodstuffs.

Sterilization

Method: Gamma irradiation Minimum Dose: 25kGy Sterility Assurance Level: 10⁻⁶

Shelf Life & Storage

Five (5) years from date of manufacture. Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.

Physical Properties

| CHARACTERISTICS | VALUE | TEST METHOD |
|-----------------------|-----------------------|---------------|
| Freedom From Holes | AQL 0.65 | EN 374-2:2003 |
| | Performance Level 3 | |
| Glove Length | 290/11.4 | EN 420:2003 |
| +/-10mm (mm/inches) | | |
| Single Wall Thickness | Cuff: 2.76/0.07 | EN 455-2:2009 |
| Minimum (mils/mm) | Palm: 3.94/0.10 | |
| | Finger Tip: 4.72/0.12 | |
| Force At Break | ≥9 N | EN 455-2:2009 |
| During Shelf Life | | |

| PALM WIDTH\$ (mm) ± 5mm | | | | | | |
|-------------------------|------|------|------|------|-----|------|
| SIZE | 5.0- | 6.0- | 7.0- | 8.0- | 9.0 | 10.0 |
| SIZL | 5.5 | 6.5 | 7.5 | 8.5 | | |
| mm | 76 | 86 | 95 | 106 | 115 | 120 |
| TEST METHOD EN 420:200 | | | | 003 | | |

Glove Cleanliness Characteristics

| | PART | ICLES | | | TYPICA | AL PARTICLE | COUNT | TEST METHOD |
|------------------------------------|----------|----------|----------|-----------|-----------------|-------------|----------|-----------------|
| ≥ 0.5µm(counts/cm²) | | | 850 | | IEST-RP-CC005.4 | | | |
| TYPICAL EXTRACTABLE IONS | | | | | TEST METHOD | | | |
| ANIONS | Fluoride | Chloride | Nitrite | Bromide | Nitrate | Phosphate | Sulphate | IEST-RP-CC005.4 |
| (µg/cm²) | ND | 0.194 | ND | ND | 0.884 | ND | 0.025 | |
| CATIONS | Lithium | Sodium | Ammonium | Potassium | Calcium | Magnesium | Zinc | |
| (µg/cm ²) | ND | 0.099 | 0.003 | 0.315 | 0.652 | 0.007 | ND | |
| ND = Not Detected, NT = Not Tested | | | | | | | | |







TO ORDER

| Re-order code | Size |
|---------------|---------|
| S-BFAP-5055 | 5.0-5.5 |
| S-BFAP-6065 | 6.0-6.5 |
| S-BFAP-7075 | 7.0-7.5 |
| S-BFAP-8085 | 8.0-8.5 |
| S-BFAP-9090 | 9.0 |
| S-BFAP-1010 | 10.0 |

Packing


One pair per inner PE wallet; one wallet per sealed EasyTear™ PE pouch; 10 pouches per sealed outer PE bag; 20 outer bags per lined carton (200 pairs)

For additional information visit us at www.bioclean.com, or call us at +44 1638 663338

Nitritex Ltd

Minton Enterprise Park Oaks Drive Newmarket, Suffolk CB8 7YY United Kingdom

S-BFAP/1806/PDS26



*Warning: No glove provides complete protection against all chemicals. Users must test the gloves against the particular chemicals and environment in which they will be used.

Ansell, ® and TM are trademarks owned by Ansell Limited or one of its affiliates. © 2012 – 2018 Ansell Limited. All Rights Reserved.



^{**}Please see product validation pack or contact Ansell customer service for specific data on use of gloves with cytotoxic drugs. Gloves used for protection against such drugs must be selected specifically for the type of chemicals used.